

Innorack® Mouse 3.5 double and single-sided.

Innorack® Mouse 3.5, dimensions and cage positions.



Innorack® Mouse 3.5, dimensions and cage positions.					
COLUMNS		6	7	8	
		CAGE POSITIONS¹			HEIGHT²
ROWS	9	108	126	144	64.25"
	10	120	140	160	70.50"
	11	132	154	176	76.75"
	12	144	168	192	83"
LENGTH		62.83"	72.13"	81.43"	

Mouse rack system includes: Innorack IVC Mouse, trolley, supply and exhaust blowers, digital controller, and power supply.

* Depth: 31.3" for double-sided; 18" for single-sided rack.

¹ Single-sided rack has half the number of cage positions shown above.

² Add 12.5" to height for blowers.



Innorack® Rat 3.5, double and single-sided.

Innorack® Rat 3.5, dimensions and cage positions.



Innorack® Rat 3.5, dimensions and cage positions.					
COLUMNS		3	4	5	
		CAGE POSITIONS ¹			HEIGHT ²
ROWS	6	N/A	48	60	60"
	7	N/A	56	70	68.75"
	8	48	64	80	77.5"
LENGTH		49.06"	62.66"	76.25"	

Rat rack system includes: Innorack IVC Rat, trolley, supply and exhaust blowers, digital controller, and power supply.

* Depth: 36.3" for double-sided; 19.8" for single-sided rack.

¹ Single-sided rack has half the number of cage positions shown above.

² Add 14.75" to height for blowers.



Innorack® IVC Tall Rat, double and single-sided.

Innorack® IVC Tall Rat, dimensions and cage positions.



Innorack® IVC Tall Rat, dimensions and cage positions.				
COLUMNS		4	5	
		CAGE POSITIONS ¹		HEIGHT ²
ROWS	5	N/A	50	63.50"
	6	48	60	74.50"
LENGTH		62.66"	76.25"	

Rat rack system includes: Innorack IVC Tall Rat, trolley, supply and exhaust blowers, digital controller, and power supply.

* Depth: 36.3" for double-sided; 19.8" for single-sided rack.

¹ Single-sided rack has half the number of cage positions shown above.

² Add 14.75" to height for blowers.



Innorack® IVC Tall Rat XL, double and single-sided.

Innorack® IVC Tall Rat XL, dimensions and cage positions.

Innorack® IVC Tall Rat XL, dimensions and cage positions.				
COLUMNS		1	2	
		CAGE POSITIONS¹		HEIGHT²
ROWS	6	12	24	74.50*
	LENGTH	62.66*	62.66*	

Rat rack system includes: Innorack IVC Tall Rat XL , trolley, supply and exhaust blowers, digital controller, and power supply.

* Depth: 36.3" for double-sided; 19.8" for single-sided rack.

¹ Single-sided rack has half the number of cage positions shown above.

² Add 14.75" to height for blowers.